

Technical sheet

Idylla™ MSI Assay

RUO

The **Idylla™ MSI Assay**, is intended for the qualitative detection of a novel panel of seven monomorphic homopolymer biomarkers for identification of human cancers with microsatellite instability (MSI). The **Idylla™ MSI Assay** uses formalin-fixed, paraffin-embedded (FFPE) tissue sections from human cancer tissue, from which nucleic acids are liberated, PCR amplified and then analyzed by high-resolution melting detection. The **Idylla™ MSI Assay** automates the entire process from FFPE sample preparation to reporting of MSI status.



Features

Idylla™ MSI Biomarkers	
ACVR2A	2q22.3-q23.1
BTBD7	14q32.12
DID01	20q13.33
MRE11	11q21
RYR3	15q13.3-q14
SEC31A	4q21.22
SULF2	20q13.12
Sample Processing Control	The MSI specific software or TTP will automatically check the validity of the measured fluorescence profiles: presence of specific PCR amplicons will result in biomarker-specific fluorescence profiles, which eliminates the need for an additional sample processing control in the cartridge.
Specimen requirements	
Sample Type	FFPE tissue sections (5-10 µm)
Neoplastic cells	≥20%, if less, macro dissection is required
Tissue area	50-600mm² (5 µm) 25-300mm² (10 µm)
Total turnaround time	
Time	150 minutes

Idylla™ MSI Posters & Abstracts

- F De Craene B. et al. Idylla™ MSI in gastric samples. ESMO 2017 poster 697P
- F Maertens G. et al. Idylla™ MSI in CRC. ESMO 2017 poster 138P
- F De Craene B. et al. Detection of microsatellite instability (MSI) in colorectal cancer samples with a novel set of highly sensitive markers by means of the Idylla™ MSI Test prototype. Online ASCO Abstract #e15639
- F Zhao et al. Association of a novel set of 7 homopolymer indels for detection of MSI with tumor mutation burden and total indel load in endometrial and colorectal cancers. Online ASCO Abstract #e15654

Disclaimer

Biocartis and Idylla™ are registered trademarks in Europe, the United States and other countries. The Biocartis trademark and logo and the Idylla trademark and logo are used trademarks owned by Biocartis. Idylla™ platform is CE-marked IVD in Europe. Idylla™ MSI Assay is available for Research Use Only (RUO), not for use in diagnostic procedures. Idylla™ is available for sale in EU, USA and some other countries. Please check availability with the local Biocartis sales representative.